



University of California  
San Francisco

# Multiplexed Amplicons for Drugs, Diagnostics, Diversity and Differentiation using High Throughput Targeted Resequencing (MAD<sup>4</sup>HaTTeR)

## Standard Operating Procedure



Adapted from CleanPlex amplicon sequencing technology by Paragon Genomics.  
Version 2 - January 2025

## **Primer pools:**

MAD<sup>4</sup>HaTTeR is a modular assay. Samples can be processed with either or both of two modules: Diversity or Resistance<sup>+</sup>.

### **Diversity:**

The diversity module consists of loci of high heterozygosity in available whole genome data.

The diversity module can be amplified with primer pool D1.1 (previously 1A)\*

### **Resistance<sup>+</sup>:**

The resistance<sup>+</sup> module consists of loci in drug resistance markers, the *hrp2/3* genes associated with diagnostic resistance, and other loci of interest for immune response studies (e.g. *csp*).

The full resistance<sup>+</sup> module should be amplified with 2 primer pools (R1.1 and R2.1) that are complementary and incompatible in a multiplexed PCR reaction due to loci overlapping. Pool R1.2 is a subset of R1.1 that significantly reduces primer dimer formation and is also incompatible with pool R2.1. Pools R1.1 and R1.2 contain most of the molecular markers of interest, and thus should take precedence over pool R2.1 if only one mPCR reaction is desired.

The EPPIcenter currently recommends using primer pools R1.2 and R2.1 for drug and diagnostic resistance, and *csp*.

A summary of the currently tested primer pools can be found here:

Module	Pool	Intended reaction	Version	Previous name	Incompatible with (overlaps)	Subset of	Number of targets	Contents
Diversity	D1.1	1	1	1A			170	Targets with high heterozygosity and Plasmodium spp. identification
Resistance <sup>+</sup>	R1.1	1	1	1B	R2.1, R1.2		82	Drug and diagnostic resistance, Plasmodium spp. identification, immune-related targets
Resistance <sup>+</sup>	R2.1	2	1	2	R1.1, R1.2		31	Drug and diagnostic resistance, Plasmodium spp. identification, immune-related targets. Complements R1.1 for missing codons
Resistance <sup>+</sup>	R1.2	1	2	5	R2.1, R1.1	R1.1	47	Drug and diagnostic resistance, Plasmodium spp. identification, immune-related targets. Complements R1.1 for missing codons, subset of R1.1 to reduce primer dimer and increase sensitivity

A full list with primer sequences, genomic locations, and genes can be found in [eppicenter.ucsf.edu/resources](http://eppicenter.ucsf.edu/resources) in the “Pool details” link in the MAD<sup>4</sup>HaTTeR section.

## **Naming convention**

As of January 2025, pools have been renamed to allow for consistency and future versions. The current nomenclature will be followed: The main use case is denoted with a letter, the intended reaction (for incompatible pools) with a number, followed by versions of those pools indicated by another number after a period. The corresponding old and new nomenclature can be found in the table above.

Diversity pools 3 and 4 were subsets of D1.1. Pool 3: 153 targets with high heterozygosity and Plasmodium spp. identification. Pool 4: 50 targets with high heterozygosity and Plasmodium spp. Identification, a subset of D1.1 to increase throughput. These pools have never been used and thus will not be assigned names under the new convention

## **PCR programs**

Save the following programs in a thermal cycler before starting

### **mPCR (multiplex PCR)**

- Initial denaturation: 95 °C 00:10:00
- Amplification for X number of cycles (**refer to table below**):
  - Denaturation: 98 °C 00:00:15 with ramping 3 °C/s
  - Annealing/Extension: 60 °C 00:05:00 with ramping 2 °C/s
- Hold at 10 °C
- Volume: 10 µL
- Lid: 105 °C

<b>Parasitemia</b>	<b>Total number of mPCR cycles</b>
≥100 p/µL	15
<100 p/µL	20

We recommend saving 2 programs with different names (mPCR15 and mPCR20) and protecting them with password or read-only settings

### **Digestion**

- 37 °C infinite hold\*
- 37 °C for 00:10:00
- Volume: 20 µL
- Lid: OFF

\* The infinite hold is used to pre-heat the block to 37 °C to have it ready to start the digestion when the samples are mixed with the digestion reagent. In some thermal cyclers, the infinite hold can be exited with a “next step” function. You will need to figure out how your instrument works.

We recommend saving the program and protecting it with password or read-only settings

### **iPCR (indexing PCR)**

- Initial denaturation: 95 °C 00:10:00
- Amplification for 15 cycles:
  - Denaturation: 98 °C 00:00:15 with ramping 3 °C/s
  - Annealing/Extension: 60 °C 00:01:25 with ramping 2 °C/s
- Hold at 10 °C
- Volume: 40 µL
- Lid: 105 °C

We recommend saving the program and protecting it with password or read-only settings

## **Personal protective equipment**

PPE should be used in accordance with your laboratory regulations. Proper use of PPE also minimizes sample contamination. Special attention should be paid to glove usage and we recommend changing them frequently

## **Multiplexed PCR (mPCR)**

### **Objective**

Amplify target regions from *P. falciparum* genomes in 1 or 2 reactions. 2 reactions are needed if using incompatible primer pools, and they can be combined once the targets have been amplified

### ***In Laminar Flow Hood in Pre-PCR room***

- Thaw primer pools at room temperature. Keep on ice after thawing
- If expecting to run full protocol in one day, thaw CP Reagent Buffer (white tube) and keep on ice after thawing
- Bring 5X mPCR Master Mix (green tube) into PCR Workstation in a cold rack
- Combine the following volumes to prepare the **mPCR reaction mix**. Keep mix on ice or cold rack. Vortex reagents to mix and briefly spin down before opening.

**(Generic recipe can be found in the next page)**

**The following recipes are for 2 reactions/sample, one with pools D1.1+R1.2 (or D1.1+R1.1, ...) and one with R2.1.**

### **For primer pools D1.1+R1.2**

<b>Reagent</b>	<b>Volume* (for <math>\geq 100</math> p/<math>\mu</math>L)</b>	<b>... x # reactions x 1.10 (10% extra**)</b>
5X mPCR Master Mix (green tube)	2 $\mu$ L	
Primer pool D1.1	0.5 $\mu$ L	
Primer pool R1.2	0.5 $\mu$ L	
Nuclease-free H <sub>2</sub> O If you are using less than 6 $\mu$ L of input DNA, increase water volume (e.g. if using 3 $\mu$ L DNA, add water up to 7 $\mu$ L)	1 $\mu$ L (or 4 $\mu$ L if using 3 $\mu$ L input DNA)	

\*: For **<100 p/ $\mu$ L** use 0.25  $\mu$ L of each primer. The total volume of water would be 1.5 or 4.5  $\mu$ L

\*\* : Depending on your setup and expected dead volumes you may want to factor up to 15% extra

### **For primer pool R2.1**

<b>Reagent</b>	<b>Volume* (for <math>\geq 100</math> p/<math>\mu</math>L)</b>	<b>... x # reactions x 1.10 (10% extra**)</b>
5X mPCR Master Mix (green tube)	2 $\mu$ L	
Primer pool R2.1	0.5 $\mu$ L	
Nuclease-free H <sub>2</sub> O If you are using less than 6 $\mu$ L of input DNA, increase water volume (e.g. if using 3 $\mu$ L DNA, add water up to 7 $\mu$ L)	1.5 $\mu$ L (or 4.5 $\mu$ L if using 3 $\mu$ L input DNA)	

\*: For **<100 p/ $\mu$ L** use 0.25  $\mu$ L of each primer. The total volume of water would be 1.5 or 4.5  $\mu$ L.

\*\* : Depending on your setup and expected dead volumes you may want to factor up to 15% extra

**Generic recipe (add 10% for dead volumes):**

Reagent	Volume * (for $\geq 100$ p/ $\mu$ L)	... x # reactions x 1.10 (10% extra)
5X mPCR Master Mix (green tube)	2 $\mu$ L	
Each primer pool	0.5 $\mu$ L	
Nuclease-free H <sub>2</sub> O If you are using less than 6 $\mu$ L of input DNA, increase water volume (e.g. if using 3 $\mu$ L DNA, add water up to 7 $\mu$ L)	Up to 4 $\mu$ L	

\*: For **<100 p/ $\mu$ L** use 0.25  $\mu$ L of each primer. The total volume of water would be 1.5 or 4.5  $\mu$ L

\*\* : Depending on your setup and expected dead volumes you may want to factor up to 15% extra

**Note: pools R1.1/R1.2 and R2.1 are incompatible.**

- Vortex mixes, briefly spin and keep on ice
- Aliquot 4  $\mu$ L (or 7  $\mu$ L if using 3  $\mu$ L input DNA) mPCR mix into PCR tubes/wells (single tubes, strips or plate).<sup>1</sup>  
Keep tubes or plate on ice or cold rack
- Put primers and 5X mPCR Master Mix back in freezer
- Transfer tubes to **Laminar Flow Hood in Pre-PCR Room**

***In Laminar Flow Hood in Pre-PCR room***

- Add 6  $\mu$ L (or 3  $\mu$ L) DNA sample to each labeled tube/well, independent of parasitemia.  
Note: total volume is 10  $\mu$ L  
Keep tubes or plate on ice or cold rack
- Vortex and briefly spin down to collect liquids and keep on ice or cold rack
- Transfer the tubes to a thermal cycler in Post-PCR room. Ice buckets or cold racks should not enter the Post-PCR room.
- Double-check that the settings are correct and run the corresponding pre-set **mPCR** program:  
Initial denaturation: 95 °C 00:10:00  
Amplification for X number of cycles (**15 for  $\geq 100$  p/ $\mu$ L, 20 for  $< 100$  p/ $\mu$ L**):  
    Denaturation: 98 °C 00:00:15 with ramping 3 °C/s  
    Annealing/Extension: 60 °C 00:05:00 with ramping 2 °C/s
  - Hold at 10 °C
  - Volume: 10  $\mu$ L
  - Lid: 105 °C

<sup>1</sup> You may do this 1 tube/well at a time, or transfer the volume into a reservoir or PCR tube strip and use a multichannel

## Preparation of reaction mixes

### After starting the thermocycler protocol prepare for next steps

- Bring Clean Mag Magnetic Beads and TE buffer to room temperature
- Bring STOP buffer to room temperature and aliquot into PCR tube strip (~200 µL per tube) so that you can use a multichannel
- Make 80% ethanol with nuclease-free water (you will need 800 µL per sample, or 400 µL if stopping after the first bead clean-up)

### If you are not stopping at the safe stopping point after the first bead clean-up, also do this:

- Bring index primers out of the freezer and thaw on ice
- Make a plan for sample indexing. Write down what index you will use for each sample in the table at the end of this worksheet
- Make digestion and indexing PCR mixes

#### *In DNA-free Laminar Flow Hood in Pre-PCR room*

If you are splitting the protocol in 2 days, make these 2 mixes at the beginning of the second day

- Combine the following volumes to prepare the **Digestion reaction mix:**  
Vortex reagents to mix and briefly spin down before opening.

Reagent	Volume	... x # reactions x 1.10 (10% extra*)
Nuclease free water	6 µL	
CP reagent buffer	2 µL	
CP digestion reagent	2 µL	
(optional) TE buffer**	10 µL	

\*: Depending on your setup and expected dead volumes you may want to factor up to 15% extra

\*\* : If you are stopping after the safe stopping point, do not add TE buffer.

**Note: Make sure that CP Reagent Buffer is completely dissolved before using**

- Combine the following volumes to prepare the **indexing PCR reaction mix:**  
Vortex reagents to mix and briefly spin down before opening.

Reagent	Volume	... x # reactions x 1.10 (10% extra*)
Nuclease free water	18 µL	
5X 2 <sup>nd</sup> PCR Master Mix	8 µL	
TE buffer**	10 µL	

\*: Depending on your setup and expected dead volumes you may want to factor up to 15% extra

\*\* : We recommend mixing TE buffer into the reaction mix and resuspending beads directly into this mixture. The original protocol called for resuspension in TE and then adding the reaction mix.

- Put CP reagent buffer and digestion reagent, and 2<sup>nd</sup> PCR Master Mix back in freezer

## Multiplexed PCR (mPCR), continued

### *In bench in Post-PCR room*

- Once mPCR reaction is done, remove tubes from the thermal cycler
- **Immediately** add STOP buffer and mix reactions if needed:

**For 2 reactions** (e.g. D1.1+R1.2 and R2.1):

Mix the 2 reactions for each sample. **Be careful not to mix different samples together**

- Add 4  $\mu\text{L}$  STOP buffer to one of the reactions for each sample
- Change the volume in the pipet to 14  $\mu\text{L}$  (or 10  $\mu\text{L}$  if that's the maximum)<sup>2</sup> and transfer all the volume to the corresponding sample tubes for the other reaction so that both reactions and the STOP buffer are combined in one tube.

Note: total volume is 24  $\mu\text{L}$

**For 1 reaction** (e.g. D1.1.+R1.1 only)<sup>3</sup>:

- Add 2  $\mu\text{L}$  STOP buffer to each of the tubes
- Add 10  $\mu\text{L}$  TE Buffer to each of the tubes

Note: total volume is 22  $\mu\text{L}$

## Post-mPCR Bead Clean-up

### **Objective**

Remove mPCR reaction leftovers (enzymes, salts and unused primers) and, more importantly, primer dimers from the solution

### *In bench in Post-PCR*

- If you haven't, prepare 80% ethanol with nuclease-free water
- Make sure that CleanMag Magnetic Beads are at room temperature and are well mixed
- Add 1.3X bead solution to each tube:
  - 31  $\mu\text{L}$  if you mixed 2 reactions
  - 29  $\mu\text{L}$  if you only had 1 reaction per sample
- Vortex vigorously to mix, or pipet up and down, and spin to collect liquids (not long enough to pellet beads)

**After this step, and until resuspension in TE Buffer, do not vortex the mixture and treat it carefully**

- Incubate for 00:05:00 minutes at room temperature
- Place on magnetic stand for 00:03:00 minutes or until the beads are collected on the side of the tubes/wells and the liquid is clear
- Remove all the liquid with a pipet set to >60  $\mu\text{L}$ . Here and in all the following steps, avoid scraping the tube wall where the beads are.
- Briefly spin down\*, place back in the magnetic stand and remove the liquid leftovers using a P20 or P10 pipet set to maximum volume. Choose a pipet with thin pipet tips to better take leftovers.

**\*: Place the tubes in the spinner with the beads so that they are on the outside, further away from the center or axis of rotation so that centrifuge force doesn't push them towards the opposite wall**

<sup>2</sup> You may use the same tips to add STOP buffer and move and combine the 2 reactions for the same sample

<sup>3</sup> To use less tips, you may combine TE and STOP buffers and add together

- Add 180  $\mu$ L 80 % ethanol
- To wash the beads, rotate the tubes so that the beads migrate from one wall to the other. In magnetic stand for plates with vertical magnets, move the plate or strip by one column.
- Incubate for 00:02:00 or until all beads have migrated. You may need to **CAREFULLY** flick the tubes/wells
- Remove the liquid with a P200 pipet
- Repeat wash: add 180  $\mu$ L 80% ethanol
- Rotate or move the tubes/plate so that the beads migrate from one wall to the other.
- Incubate for 00:02:00 or until all beads have migrated. You may need to **CAREFULLY** flick the tubes/wells
- Remove the liquid with a P200 pipet
- Briefly spin down, place back in the magnetic stand and remove all the remaining liquid using a P20 or P10 pipet (this is only on the second wash, not the first one)
- Leave tubes/wells open to dry at room temperature. **Generally, a 5 min dry time is enough, but it will depend on room temperature and humidity.** You may not have to wait and just add TE buffer as soon as you finish removing the liquid from the last well  
The beads should look matte, not shiny. **Under-drying (carrying ethanol) and over-drying (cracking) can lead to reduced yield**
- Add 10  $\mu$ L TE buffer (or 20  $\mu$ L of the mixture TE buffer and digestion reaction mix if proceeding to the next step immediately) equilibrated at room temperature. Close the tubes/wells and vortex vigorously and/or tap to re-suspend the beads.<sup>4</sup>
- Quickly spin down to collect liquids  
The magnetic beads will not affect the rest of the reactions, there is no need to remove them

~~~~~SAFE STOPPING POINT~~~~~

**If you want to stop here, store at -20 °C**

---

<sup>4</sup> If you are continuing with the full protocol resuspend directly in a mixture of TE and the digestion reaction mix



**If you stopped in the previous step and left samples at -20 °C, at the beginning of day 2:**

- Thaw CP Reagent Buffer (white tube) and keep on ice after thawing
- Bring index primers out of the freezer
- Bring Stop Buffer, TE Buffer and CleanMag Magnetic Beads to room temperature
- Make 80% ethanol with nuclease-free water (400 µL per sample)
- Write down what index you will use for each sample in the table at the end of this worksheet
- Make digestion and indexing PCR mixes in DNA-free Laminar Flow Hood
- Bring samples out of the freezer, thaw and spin down to collect liquids

## **Digestion**

### **Objective**

Remove non-specific products that are not double-stranded DNA

Warning: if the digestion reaction is not stopped in time, degradation of library products can occur. Promptly add STOP buffer and proceed to bead clean-up.

#### ***In bench in Post-PCR room***

- Start the protocol without the samples. Once the block has reached 37 °C press Pause
- Vortex Digestion reaction mix and briefly spin down to collect liquids
- Add 10 µL of digestion reaction mix to each tube<sup>5</sup>
- Vortex and quickly spin down to collect liquids
- Transfer tubes to thermal cycler pre-set at 37 °C and press Resume (or equivalent option in your thermal cycler) to incubate for 10 min
- After 10 min, **immediately** add 2 µL of Stop Buffer in each tube

## **Post-Digestion Bead Clean-up**

### **Objective**

Remove enzymes, salts, and short digestion products from the solution

#### ***In bench in Post-PCR room***

- If you haven't, prepare 80% ethanol with nuclease-free water
- Make sure that CleanMag Magnetic Beads are at room temperature
- Add 1.3X bead solution to each tube: 29 µL
- Vortex vigorously, or pipet up and down, to mix and spin to collect liquids (not long enough to pellet beads)

**After this step, and until resuspension in TE Buffer, do not vortex the mixture and treat it carefully**

- Incubate for 00:05:00 minutes at room temperature
- Place on magnetic stand for 00:03:00 minutes or until the beads are collected on the side of the tubes/wells and the liquid is clear
- Remove all the liquid: first with a pipet set to >60 µL
- Briefly spin down\*, place back in the magnetic stand and remove the liquid leftovers using a P20 or P10 pipet set to maximum volume. Choose a pipet with thin pipet tips to better take leftovers.

**\*: Place the tubes in the spinner with the beads so that they are on the outside, further away**

---

<sup>5</sup> If you did not stop at the previous safe stopping point, add the digestion reaction mix together with TE buffer in the last steps of the previous section

**from the center or axis of rotation so that centrifuge force doesn't push them towards the opposite wall**

- Add 180 µL 80 % ethanol
- To wash the beads, rotate the tubes so that the beads migrate from one wall to the other. On magnetic stand for plates with vertical magnets, move the plate or strip by one column.
- Incubate for 00:02:00 or until all beads have migrated. You may need to **CAREFULLY** flick the tubes/wells
- Remove the liquid with a P200 pipet
- Repeat wash: add 180 µL 80% ethanol
- Rotate or move the tubes/plate so that the beads migrate from one wall to the other.
- Incubate for 00:02:00 or until all beads have migrated. You may need to **CAREFULLY** flick the tubes/wells
- Remove the liquid with a P200 pipet
- Briefly spin down, place back in the magnetic stand and remove all the remaining liquid using a P20 or P10 pipet (this is only on the second wash, not the first one)
- Leave tubes/wells open to dry at room temperature. **Generally, a 5 min dry time is enough, but it will depend on room temperature and humidity.** You may not have to wait and just add the mixture of TE and indexing PCR mix as soon as you finish removing the liquid from the last well  
The beads should look matte, not shiny. **Under-drying (carrying ethanol) and over-drying (cracking) can lead to reduced yield**
- Resuspend directly into 36 µL of mixture of TE and indexing PCR mix, 36 µL Close the tubes/wells and vortex vigorously and/or tap to re-suspend the beads.
- Quickly spin down to collect liquids

## Indexing PCR

### Objective

Further amplify amplicons and append adaptors so that the DNA binds to the sequencer flow cell, and a unique DNA barcode in each sample

#### *In bench in Post-PCR room*

- Vortex Indexing PCR mix and briefly spin down to collect liquids
- Add 26 µL Indexing PCR mix to each tube
- Spin down indexing primer plate and cut out the cover from the wells you will use
- Add 4 µL of indexing primers from a unique index plate well into each of the samples
- Reseal indexing PCR plate
- Vortex and briefly spin to collect liquids
- Double-check that the settings are correct and run the pre-set **iPCR** program:  
Initial denaturation: 95 °C 00:10:00  
Amplification for 15 cycles :  
    Denaturation: 98 °C 00:00:15 with ramping 3 °C/s  
    Annealing/Extension: 60 °C 00:01:15 with ramping 2 °C/s
  - Hold at 10 °C
  - Volume: 40 µL
  - Lid: 105 °C
- Make note of the primers used
- After the program is done, take tubes out and store or keep using

~~~~~SAFE STOPPING POINT (store at -20°C)~~~~~

## **Capillary electrophoresis**

### **Objective**

Assess the quality of libraries:

1. Do we have indexed amplicons?
2. Do we have primer dimers?
3. What are our library concentrations?

- Briefly spin down tubes and place on magnetic stand to separate beads, which cannot be loaded into capillary electrophoresis systems
- Select random samples to run on capillary electrophoresis. Include negative and positive controls when available
- Follow the instructions for the specific instrument you are using

### **Results interpretation:**

1. Internal controls: All samples should have lower and upper markers that come with the buffer. These markers are used to back-calculate the length of DNA fragments based on the time they were detected. If no markers are observed, the readings should not be interpreted
2. Negative control: A negative control with no template for mPCR should still have primer dimers that are amplified and are indexed. A peak should be observed at ~200 bp and no peak should be observed at ~400 bp
3. Positive control (optional) and samples: Positive control and samples should have an amplified library (~400 bp). All of these samples will also have a peak at ~200 bp corresponding to the primer dimers that were not eliminated during bead cleanups and were amplified in indexing PCR. Library concentrations would optimally be at ~10 nM.

## **Sample pooling**

### **Objective**

Combine libraries from multiple samples to be sequenced

### ***In bench in Post-PCR room***

- Before starting make 80% ethanol with nuclease-free water and bring CleanMag Magnetic Beads to room temperature
- Create a sample sheet (you may use the templates available at [epicenter.ucsf.edu/resources](http://epicenter.ucsf.edu/resources))
- Double check that indexes and sample names are not duplicated
- Briefly spin down tubes and place on magnetic stand to separate beads
- Pool samples by mixing them into a single 1.5 mL microcentrifuge tube

We recommend using the following volumes if using 15 cycles in mPCR:

- 30  $\mu$ L for 1 p/ $\mu$ L
- 20  $\mu$ L for 10 p/ $\mu$ L
- 15  $\mu$ L for 100 p/ $\mu$ L
- 6  $\mu$ L for 1,000 p/ $\mu$ L
- 3  $\mu$ L for 10,000 p/ $\mu$ L

If you have samples with unknown parasitemia, but suspect they are high, mix 10  $\mu$ L of each sample.

If you have capillary electrophoresis for each of the samples, pool with volumes inversely proportional to the concentration of the 300-500 bp region

**Total volume in pool: .....**

## Post-Pooling Bead cleanup

### Objective

Remove leftover primer dimers and indexing PCR reagents

- Add 1X volume of beads to the pool (same volume you entered above). You may need to split in more than 1 tube if the total volume is > 1.5 mL
- Vortex and briefly spin to collect liquids (but don't pellet beads)
- Incubate for 00:05:00 at room temperature
- Place tube(s) in microcentrifuge tube magnetic stand and incubate for 00:03:00 or until the solution is clear
- Remove supernatant
- Briefly spin down and place back in magnetic stand, remove the rest of the supernatant
- Add 1.5 mL 80% ethanol
- Rotate tube(s) and incubate until all beads have migrated to the opposite wall
- Remove ethanol supernatant
- Add 1.5 mL 80% ethanol
- Rotate tube(s) and incubate until all beads have migrated to the opposite wall
- Remove ethanol supernatant
- Briefly spin down, place back in magnetic stand and remove the rest of the supernatant
- Leave tube(s) open to dry at room temperature. **Generally, a 5 min dry time is enough, but it will depend on room temperature and humidity.** The beads should look matte, not shiny. **Under-drying (carrying ethanol) and over-drying (cracking) can lead to reduced yield**
- Resuspend beads with 43  $\mu$ L TE Buffer.  
If using multiple tubes, resuspend in one tube and use that resuspension to resuspend the rest of the tubes so that the final total volume is 43  $\mu$ L
- Incubate at Room Temperature for 00:02:00
- Place tube in magnetic stand and incubate at room temperature for 00:03:00 or until liquid is clear
- Transfer 40  $\mu$ L to a clean, labeled tube

## Capillary electrophoresis

### Objective

Assess pool quality and determine need of further purification

- Follow the instructions for the specific instrument you are using

### Results interpretation and recommendations:

The following concentrations, volumes and percentages are points of reference, and could be modified depending on laboratory needs (sequencing requirements, general yield, etc). At the EPPIcenter we have better yields with gel purification than bead cleanup so our preferred method for high dimer contents is gel purification. If the library concentration is high enough, we prefer a bead cleanup for simplicity.

| <b>Primer dimer</b> | <b>400 bp concentration</b> | <b>Recommendation</b>   |
|---------------------|-----------------------------|---|
|                     |                             | For all, assess quality with capillary electrophoresis and follow these recommendations again   |
| >5%                 | >20 nM                      | Perform a 1X bead clean up and resuspend in 15 µL   |
| >5%                 | <20 nM                      | Perform gel purification to avoid losing library and more stringently remove primer dimers  |
| 1-5%                | >50 nM                      | Perform a 1X bead clean up and resuspend in 15 µL   |
| 1-5%                | 20-50 nM                    | Perform gel purification to avoid losing library and more stringently remove primer dimers  |
| 1-5%                | <20 nM                      | Proceed with the pool as is, no further purification is recommended to avoid losing the library. You can consider performing a 0.8X bead cleanup and resuspending in 15 µL to concentrate the library. Note that most sequencing cores will require at least 4 nM |
| <1%                 | <20 nM                      | Proceed with pool as is, no further purification is recommended to avoid losing library   |
| <1%                 | >20 nM                      | Perform a 1X bead clean up and resuspend in 15 µL   |

## **Gel purification**

### **Objective**

Remove any remaining primer dimers

### **Procedure**

1. Cast a 2.5% agarose gel in TAE buffer with 1X SYBRsafe
2. Place the gel in an electrophoresis system and fill up with TAE
3. Load 5-10  $\mu\text{L}$  DNA ladder in one well
4. Add 4  $\mu\text{L}$  6X loading buffer to 20  $\mu\text{L}$  of the pooled libraries keeping half in case you need to repeat. In case the concentration of the libraries is too low (as a rough point of reference use 10 nM), consider running the full volume. Vortex and spin to collect liquids
5. Load the pool into 1 or more lanes (depending on comb size you may not be able to fit in one lane). Leave an unused lane between the ladder and the pool
6. Run with constant voltage at 140 V for 1 h
7. After 1 h, quickly image the gel. If the bands (200 and 400 bp) are clearly separated, you may continue to excise. Otherwise, run for longer
8. Excise the 400 bp band
9. Using a DNA gel extraction kit, dissolve the excised gel and run through a column following manufacturer's instructions
10. Elute with 15  $\mu\text{L}$  elution buffer

## **Capillary electrophoresis**

### **Objective**

Ensure that primer dimers are at less than 5% of the total sample

- Follow the instructions for the specific instrument you are using

### **Results interpretation:**

For the pool to be ready for sequencing, we need to check 2 measurements:

1. The primer dimer peak, if present, should not be more than 5% of the sample. If not, you will need to do another gel purification
2. The libraries should be at  $\geq 4$  nM (as required per Illumina protocols)

**Was another gel purification necessary? .....**

**Appendix 1: Summary of all sample and indexes used**

| Tube No. | Sample id | Primer pools used | mPCR cycle no. | Index used (Set and well) | Volume pooled |
|----------|-----------|-------------------|----------------|---------------------------|---------------|
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |
|          |           |                   |                |                           |               |